

P.A.P.

Day I

Prepare Slides:

- Cut sections @ 3 microns on clean glass slide.
- Place in 58-60 degree C. incubator overnight.

Day II

- Place 200 ml 1/10 PBS in a 37 degree C. incubator.
- Weigh out 0.2 gm trypsin.
- Make up albuminized 1/10 PBS:

0.1 g bovine albumin per 100 ml 1/10 PBS (This solution is used to prepare all antibody dilutions and for washing slides).

Procedure

(1) Transfer slides directly from 58-60 degree C. incubator to:

1. xylene 10 min.
2. xylene 10 min.
3. xylene 10 min.
4. xylene 10 min.
5. 95% alcohol 5 min.
6. Cold Tap H₂O 5-10 min.
7. 37 degree Tap H₂O 3-5 min.

(2) Transfer slides to trypsin solution as made up previously for 6 minutes at 37 degrees C. (Discard after use).

(3) Rinse in cold tap H₂O for 5 min.

Inhibition of Endogenous Peroxidase

(4) Transfer slides to freshly prepared 200 ml of methanol + 2 ml 30% H₂O₂ for 30 minutes.

(5) Rinse in cold tap water for 5 minutes

(6) 95% alcohol 5 min.

(7) Abs. alcohol 5 min.

(8) Chloroform 5 min.

(9) Acetone 5 min.

- (10) Cold Tap H₂O 5 min.
- (11) Albuminized PBS 5 min. minimum

Primary Antibody

- (12) Prepare dilutions of primary antibody using albuminized PBS 140 - 200 μ l per slide.
- (13) Prepare slides one by one by drying outside and leaving a film of the albuminized PBS around the section (NO AIR BUBBLES). (DO NOT LET THE SECTIONS DRY!).
- (14) Cover section with diluted primary antibody (one slide at a time) for 30 minutes in a humidification chamber.
- (15) Wash each slide carefully with 1/10 PBS (albuminized) using a squeeze bottle.
- (16) Transfer to a staining dish of the albuminized PBS for 5 min.

Secondary Antibody

- (17) Prepare dilution as follows:
for 15 slides

PBS 1/10 - 2700 μ l

Normal Human Serum - 300 μ l

Shake well and discard 200 μ l of this solution.
Add 200 μ l peroxidase conjugated rabbit immunoglobulins to mouse immunoglobulins. Mix well.
- (18) Prepare slides as for primary antibodies and add secondary antibodies for 30 minutes.
- (19) Wash with albuminized PBS and return to staining dish of PBS for 5 min.

Tertiary Antibody

- (20) Prepare dilutions as for secondary antibody using peroxidase conjugated swine immunoglobulins to rabbit immunoglobulins.
- (21) Prepare slides as before and add antibody for 30 minutes.

(22) Wash well with PBS and return to staining dish of PBS for 5 min.

Revelation of Peroxidase Staining:

(23) Prepare DAB:

Dissolve 6 mg of DAB (diaminobenzidine) in 10 ml Tris-HCL (0.05M). Add 3 μ l of 30% H₂O₂ just before use. Keep in dark. (If larger quantities are needed, use several tubes ready for H₂O₂).



10 ml



10 ml



10 ml

(24) Arrange slides in humidification chamber. Put DAB on each slide for 3 to 5 minutes. (Do not attempt to reveal too many slides at one time - maximum 10 slides).

(25) Rinse with cold tap H₂O to stop reaction.

(26) Return to a staining jar of water.

Counterstain

(27) Haematoxylin for 10 to 20 seconds.

(28) Rinse in cold tap H₂O.

(29) ~~Blue~~ in STWS for 1 min.

(30) Wash in cold tap H₂O for 5 min.

(31) Dehydrate, clear and mount.